



Resprep® 12- or 24-Port SPE Manifolds

- Use with any standard male luer-end SPE cartridges.
- Inert, Teflon® sample guides reduce cross-contamination and carryover.
- Flexible sample collection rack will accommodate a variety of receiving vessels.
- Quick vacuum-release valve for better system control.
- Individual vacuum control for each cartridge—improves reproducibility.



| Description | qty. | cat.# | |
|---|------|-------|--|
| Complete Resprep SPE Manifold, 12-Port (Includes: glass basin with built-in | | | |
| vacuum regulator, polypropylene lid with 12 individual control valves, | | | |
| 12-position collection rack, 12 Teflon sample guides, and waste container) | kit | 26077 | |
| Complete Resprep SPE Manifold, 24-Port (Includes: glass basin with built-in | | | |
| vacuum regulator, polypropylene lid with 24 individual control valves, | | | |
| 24-position collection rack, and 24 Teflon sample guides) | kit | 26080 | |

Resprep® Manifold Replacement Parts

| nesprep mannola neplacement rates | | | |
|--|--------|-------|--|
| Description | qty. | cat.# | |
| Replacement Waste Container, 12-Port | ea. | 24014 | |
| Replacement Vacuum Valve and Gauge Assembly | ea. | 24008 | |
| Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Port | ea. | 25991 | |
| Collection Rack | | | |
| Collection Rack, 12-Port | ea. | 26079 | |
| Collection Rack, 24-Port | ea. | 26082 | |
| Manifold Lid | | | |
| Replacement Manifold Lid (sample guides not included), 12-Port | ea. | 26078 | |
| Replacement Manifold Lid (sample guides not included), 24-Port | ea. | 26081 | |
| Manifold Lid Replacement Parts | | | |
| Valves, 12 or 24-Port | 2-pk. | 26083 | |
| Valve Attachment | 48-pk. | 26130 | |
| Needle Attachment | 48-pk. | 26131 | |
| Sample Guides, 12 or 24-Port | 12-pk. | 26084 | |
| *Not included in 24 part manifold (cat # 26000 or 26002) | | | |

^{*}Not included in 24-port manifold (cat.# 26080 or 26082).

Oil-Free Vacuum/Pressure Station

for 12-Port Resprep® SPE Systems

Allows switching from pressure to vacuum in a matter of seconds. Quiet, oil-free unit won't contaminate the extraction system or your samples. Provides 20" Hg (68k Pa) vacuum or 18 psi (124 kPa) pressure.


| Description | qty. | cat.# |
|---|------|-------|
| Oil-Free Vacuum/Pressure Station, 115VAC, 60Hz, US | ea. | 24002 |
| Oil-Free Vacuum/Pressure Station, 230VAC, 50Hz, Europe (CE certified) | ea. | 24003 |
| Vacuum Tubing (10ft./3m, 1/4" ID) | ea. | 24016 |

Not recommended for use with 24-port manifold.

Warranty period is one-year from date of purchase. Evaluation fee is charged for repairs out of warranty.



Resprep® Manifold Replacement Parts

























Lid components prior to assembly.




Fully assembled unit shown with Glass Vacuum Chamber (cat.# 25991) and Disk Holders (cat.# 24020)



Resprep® 6-Port Disk Manifold Lid

- Low cost option for disk extraction; fits standard 3 ³/₄" x 7 ¹/₂" glass vacuum chambers.
- Doubles sample capacity—holds 6 disks, compared to standard 12-port manifolds, which only hold 3.
- Individual vacuum control for each port improves reproducibility.
- Collection plate design secures variety of receiving vessels.
- Inert Teflon® sample guides reduce cross-contamination and carryover.
- Compatible with any standard male luer end disk holder.



^{*}Vacuum chamber not included.

Resprep® Sample Delivery System

- Compatible with Resprep[®] 1, 3, 6, and 15 mL SPE cartridges and Diskcover-47 Extraction Disk Holder (cat.# 24020).
- 6 Teflon® transfer lines ($^1/_8$ " OD x $^1/_{16}$ " ID x 36" long); each contains a different embedded color stripe, for easy sample identification.
- Specified in EPA drinking water methods.
- Tested to pH of 1 to ensure no contaminants leach from system.

Use the Resprep® Sample Delivery System to transfer large volumes of low viscosity samples directly from a bottle to a solid phase extraction cartridge or to a disk on a vacuum manifold system for extraction or cleanup. Each unit consists of 6 transfer lines with a stainless steel weight on one end, and a color-coded screw fitting and polyethylene terephthalate (PET) adapter on the opposite end.

| Description | qty. | cat.# |
|--------------------------------|-------|-------|
| Resprep Sample Delivery System | 6-pk. | 26250 |
| | | |

ChromaBLOGraphy

Topical and timely insights from top chromatographers.

Visit us at blog.restek.com



